

NPR-C antibody [HMV4274] HistoMAX™

Cat. No. GTX643510

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-1:100

Note : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Target Retrieval Solution buffer.

Not tested in other applications.

Calculated MW	60 kDa. (Note)
Product Note	Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 2% BSA, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Synthetic peptide of human BRD4.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 10 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX643510 IHC-P Image

IHC-P analysis of human follicular cyst from ovary tissue using GTX643510 NPR-C antibody [HMOV4274] HistoMAX™.


A strong membranous NPR-C staining of a subset of granulosa cells.



GTX643510 IHC-P Image

IHC-P analysis of human corpus luteum from ovary tissue using GTX643510 NPR-C antibody [HMOV4274] HistoMAX™.

A distinct membranous and perinuclear cytoplasmic NPR-C staining of corpus luteum cells.



GTX643510 IHC-P Image

IHC-P analysis of human placenta (chorion) tissue using GTX643510 NPR-C antibody [HMOV4274] HistoMAX™.

A distinct membranous and cytoplasmic NPR-C staining of decidua cells.



For full product information, images and publications, please visit our [website](#).