

SH3PXD2A antibody [HL4487]

Cat. No. GTX644054

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P (cell pellet)
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P (cell pellet)	Assay dependent

Not tested in other applications.

Product Note This antibody is specific for human SH3PXD2A protein, and it does not cross-react with human HTRA1 protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide of human SH3PXD2A.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX644054 IHC-P (cell pellet) Image

SH3PXD2A antibody [HL4487] detects SH3PXD2A protein by immunohistochemical analysis.

Sample: Paraffin-embedded transfected 293T cell FFPE Cell Pellet Blocks.

Green: SH3PXD2A stained by SH3PXD2A antibody [HL4487] (GTX644054) diluted at 1:1000.

Red: DDDDK tag antibody [GT1216] (GTX636630) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTX30920).

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).