

## SSA1 antibody

Cat. No. GTX64414

|              |                |
|--------------|----------------|
| Host         | Rabbit         |
| Clonality    | Polyclonal     |
| Isotype      | IgG            |
| Applications | WB, ICC/IF, IP |
| Reactivity   | Human, Rat     |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | 1:10 - 1:100         |
| IP                 | Assay dependent      |

Not tested in other applications.

Calculated MW 54 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.05% ProClin 300                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Recombinant fusion protein containing a sequence corresponding to amino acids 80-330 of human TRIM21 (NP_003132.2).                                                                                                        |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX64414 ICC/IF Image

ICC/IF analysis of 3T3 cells using GTX64414 SSA1 antibody.

Red : Primary antibody

Blue : DAPI

Dilution : 1:100




### GTX64414 WB Image

WB analysis of rat lung tissue lysates using GTX64414 SSA1 antibody.

Dilution : 1:1000

Loading : 25µg



### GTX64414 IP Image

IP analysis of HepG2 cell lysates using GTX64414 SSA1 antibody.


Lane 1 : Input

Lane 2 : Ctrl IgG

Lane 3 : Primary antibody

IP antibody : 0.5 µg of GTX64414

Dilution : 1:500



### GTX64414 WB Image

WB analysis of 293T cell lysates and SSA1 KO 293T cell lysates using GTX64414 SSA1 antibody.

Dilution : 1:1000

Loading : 25µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).