

## DC-SIGN antibody

Cat. No. GTX64419

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB, ICC/IF |
| Reactivity   | Human      |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | Assay dependent      |

Not tested in other applications.

Calculated MW 46 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Recombinant fusion protein containing a sequence corresponding to amino acids 265-404 of human CD209 (NP_066978.1).                                                                                                        |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX64419 ICC/IF Image

ICC/IF analysis of U-937 cells using GTX64419 DC-SIGN antibody.

Red : Primary antibody

Blue: DAPI

Dilution: 1:100



### GTX64419 WB Image

WB analysis of THP-1 cell lysates using GTX64419 DC-SIGN antibody.

Loading : 25µg

Dilution:1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).