

GPR83 antibody [HL4707]

Cat. No. GTX644708

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:60000

Not tested in other applications.

Product Note This antibody was raised against human GPR83 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human GPR83. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX644708 IHC-P Image


GPR83 antibody [HL4707] detects GPR83 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

GPR83 stained by GPR83 antibody [HL4707] (GTX644708) diluted at 1:60000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTX644708 IHC-P Image


GPR83 antibody [HL4707] detects GPR83 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

GPR83 stained by GPR83 antibody [HL4707] (GTX644708) diluted at 1:60000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTX644708 IHC-P Image

GPR83 antibody [HL4707] detects GPR83 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

GPR83 stained by GPR83 antibody [HL4707] (GTX644708) diluted at 1:60000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.



For full product information, images and publications, please visit our [website](https://www.genetex.com).