

# Enterovirus 71 3C antibody [HL4973]

**Cat. No. GTX645059**

|                     |                                   |
|---------------------|-----------------------------------|
| <b>Host</b>         | Rabbit                            |
| <b>Clonality</b>    | Monoclonal                        |
| <b>Isotype</b>      | IgG                               |
| <b>Applications</b> | WB                                |
| <b>Reactivity</b>   | Coxsackievirus A6, Enterovirus 71 |

**Package**  
100 µl, 25 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000-1:10000       |

Not tested in other applications.

### Product Note

This antibody recognize 3C protein of Enterovirus 71 and Coxsackievirus A6, and it does not cross react with 3C protein of Enterovirus D68.

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | No preservatives                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | Recombinant fragment of Enterovirus 71 3C.                                                                                                                                                                                 |
| <b>Purification</b>  | Affinity purified by Protein A.                                                                                                                                                                                            |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX645059 WB Image

Non-transfected (–) and transfected (+) 293T whole cell extracts (30 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with Enterovirus 71 3C antibody [HL4973] (GTX645059) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



### GTX645059 WB Image

Non-infected (–) and infected (+) RD whole cell extracts (5 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with Enterovirus 71 3C antibody [HL4973] (GTX645059) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).