

PRRS virus Nsp1 alpha antibody [HL5078]

Cat. No. GTX645268

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	PRRS virus

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000

Not tested in other applications.

Product Note

This antibody is able to recognize PRRSV Nsp1 alpha protein and PRRSV MD006-P1, and it does not cross-react with PRRSV Nsp1 beta protein and PRRS TSYN-P1.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of PRRS virus Nsp1 alpha.(#VR-2332)
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX645268 WB Image

Non-transfected (–) and transfected (+) 293T whole cell extracts were separated by 15% SDS-PAGE, and the membrane was blotted with PRRS virus Nsp1 alpha antibody [HL5078] (GTX645268) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX645268 WB Image

Non-infected (–) and infected (+) Marc-145 whole cell extracts (15 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with PRRS virus Nsp1 alpha antibody [HL5078] (GTX645268) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).