

MCH-1R antibody [HL5082]

Cat. No. GTX645272

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human MCH-1R Intracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Synthetic peptide encompassing a sequence within the Intracellular domain of human MCH-1R. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx645272 IHC-P Image

MCH-1R antibody [HL5082] detects MCH-1R protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

MCH-1R stained by MCH-1R antibody [HL5082] (GTx645272) diluted at 1:800.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx645272 IHC-P Image

MCH-1R antibody [HL5082] detects MCH-1R protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse hippocampal CA1.

MCH-1R stained by MCH-1R antibody [HL5082] (GTx645272) diluted at 1:800.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx645272 IHC-P Image

MCH-1R antibody [HL5082] detects MCH-1R protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

MCH-1R stained by MCH-1R antibody [HL5082] (GTx645272) diluted at 1:1000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTx645272 IHC-P Image

MCH-1R antibody [HL5082] detects MCH-1R protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

MCH-1R stained by MCH-1R antibody [HL5082] (GTx645272) diluted at 1:1000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).