

VEGF Receptor 3 antibody [HL5094]

Cat. No. GTX645284

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent

Note : The FcR Blocking Reagent is recommended for use prior to the addition of the researcher's antibody to samples that contain human Fc receptor-expressing cells.

Not tested in other applications.

Observed MW (kDa)	175 kDa.
Product Note	This antibody was raised against human VEGF Receptor 3 Intracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Intracellular domain of human VEGF Receptor 3. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 11 Page 1 of 2


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX645284 WB Image

Boiled and unboiled various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with VEGF Receptor 3 antibody [HL5094] (GTX645284) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody. Corresponding RNA expression data are based on Human Protein Atlas program.


GTX645284 ICC/IF Image

VEGF Receptor 3 antibody [HL5094] detects VEGF Receptor 3 protein by immunofluorescent analysis. Sample: HEL 92.1.7 cells were fixed in 4% paraformaldehyde at RT for 15 min and treated with Fc blocker after blocking step. Green: VEGF Receptor 3 stained by VEGF Receptor 3 antibody [HL5094] (GTX645284) diluted at 1:100. Blue: Fluoroshield™ with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).