

CHCHD3 antibody

Cat. No. GTX64550

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IHC-P	1:50 - 1:100
IP	1:50 - 1:200

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human CHCHD3 (NP_060282.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).


DATA IMAGES


GTx64550 IP Image

IP analysis of mouse liver tissue lysate using GTx64550 CHCHD3 antibody.


Antibody amount : 3µg / 600µg lysate

Dilution : 1:1000


GTx64550 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx64550 CHCHD3 antibody.

Dilution : 1:100


GTx64550 IHC-P Image

IHC-P analysis of mouse lung tissue using GTx64550 CHCHD3 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).