

Galactosylceramidase antibody

Cat. No. GTX64559

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 77 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.05% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human GALC (NP_000144.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX64559 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX64559 Galactosylceramidase antibody.
Dilution : 1:100




GTX64559 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX64559 Galactosylceramidase antibody.
Dilution : 1:100



GTX64559 IHC-P Image

IHC-P analysis of human liver cancer tissue using GTX64559 Galactosylceramidase antibody.
Dilution : 1:100



GTX64559 WB Image

WB analysis of various sample lysates using GTX64559 Galactosylceramidase antibody.
Dilution : 1:1000
Loading : 25µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).