

Glycophorin A antibody [HL5297]

Cat. No. GTX645762

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	1:100-1:1000

Note : The FcR Blocking Reagent is recommended for use prior to the addition of the researcher's antibody to samples that contain human Fc receptor-expressing cells.

Not tested in other applications.

Product Note This antibody was raised against human Glycophorin A Intracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Synthetic peptide encompassing a sequence within the Intracellular domain of human Glycophorin A. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 13 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTX645762 IHC-P Image

Glycophorin A antibody [HL5297] detects Glycophorin A protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tonsil.

Glycophorin A stained by Glycophorin A antibody [HL5297] (GTX645762) diluted at 1:200.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min


GTX645762 IHC-P Image


Glycophorin A antibody [HL5297] detects Glycophorin A protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tissues.

Glycophorin A stained by Glycophorin A antibody [HL5297] (GTX645762) diluted at 1:200.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min

Corresponding RNA levels (nTPM) in the tissues are based on Human Protein Atlas program.


GTX645762 ICC/IF Image

Glycophorin A antibody [HL5297] detects Glycophorin A protein by immunofluorescent analysis.

Sample: HEL 92.1.7 cells were fixed in ice-cold MeOH for 5 min and treated with Fc blocker after blocking step.

Green: Glycophorin A stained by Glycophorin A antibody [HL5297] (GTX645762) diluted at 1:100.

Blue: Fluoroshield™ with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](#).