

Obestatin antibody [HL5397]

Cat. No. GTX645966

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000
IHC-Fr	Assay dependent

Not tested in other applications.

Product Note Highly recommended for IHC-P in mouse tissues.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide of human Obestatin.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx645966 IHC-P Image


Obestatin antibody [HL5397] detects Obestatin protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

Obestatin stained by Obestatin antibody [HL5397] (GTx645966) diluted at 1:250.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTx645966 IHC-P Image


Obestatin antibody [HL5397] detects Obestatin protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tissues.

Obestatin stained by Obestatin antibody [HL5397] (GTx645966) diluted at 1:250.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.



GTx645966 IHC-P Image

Obestatin antibody [HL5397] detects Obestatin protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat tissues.

Obestatin stained by Obestatin antibody [HL5397] (GTx645966) diluted at 1:250.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx645966 IHC-Fr Image

Obestatin antibody [HL5397] detects Obestatin protein by immunohistochemical analysis.

Sample: Frozen-sectioned mouse stomach.

Green: Obestatin stained by Obestatin antibody [HL5397] (GTx645966) diluted at 1:1000.

Red: Vimentin antibody [GT812] (GTx629744) diluted at 1:1000.

Blue: Hoechst 33342 staining.



For full product information, images and publications, please visit our [website](https://www.genetex.com).