

Avian Influenza A virus H7N7 HA (Hemagglutinin) antibody [HL5486]

Cat. No. GTX646150

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Avian Influenza A virus

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000

Not tested in other applications.

Product Note

This antibody is specific for Influenza A virus H7N7 Hemagglutinin protein (A/duck/Taiwan/33/1993), and it does not cross-react with Hemagglutinin proteins of Influenza A virus (H1N1 and H3N2), Avian Influenza A virus (H5N1, H5N3, H5N8, H9N2, and H10N3), or Influenza B virus.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide of Avian Influenza A virus H7N7 HA (Hemagglutinin) (A/duck/Taiwan/33/1993(H7N7))
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

