

Retinoid X Receptor alpha antibody [HL5604]

Cat. No. GTX646286

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	1:100-1:1000

Not tested in other applications.

Observed MW (kDa) 60 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide of human Retinoid X Receptor alpha.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX646286 IHC-P Image

Retinoid X Receptor alpha antibody [HL5604] detects Retinoid X Receptor alpha protein by immunohistochemical analysis.

Sample: Paraffin-embedded human pancreas.

Retinoid X Receptor alpha stained by Retinoid X Receptor alpha antibody [HL5604] (GTX646286) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX646286 IHC-P Image


Retinoid X Receptor alpha antibody [HL5604] detects Retinoid X Receptor alpha protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tissues.

Retinoid X Receptor alpha stained by Retinoid X Receptor alpha antibody [HL5604] (GTX646286) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.



GTX646286 IHC-P Image

Retinoid X Receptor alpha antibody [HL5604] detects Retinoid X Receptor alpha protein by immunohistochemical analysis.

Sample: Paraffin-embedded human kidney.

Retinoid X Receptor alpha stained by Retinoid X Receptor alpha antibody [HL5604] (GTX646286) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX646286 IHC-P Image

Retinoid X Receptor alpha antibody [HL5604] detects Retinoid X Receptor alpha protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

Retinoid X Receptor alpha stained by Retinoid X Receptor alpha antibody [HL5604] (GTX646286) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).