

FKBP4 antibody

Cat. No. GTX64678

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:100
IHC-P	1:50 - 1:200
IP	1:50 - 1:200

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-459 of human FKBP4 (NP_002005.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX64678 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX64678 FKBP4 antibody.
Blue : DAPI


GTX64678 WB Image

WB analysis of various sample lysates using GTX64678 FKBP4 antibody.
Dilution : 1:1000
Loading : 25µg per lane


GTX64678 IP Image

IP analysis of 293T cell lysate using GTX64678 FKBP4 antibody.
Antibody amount : 1µg / 200µg lysate
Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).