

DRIL1 antibody

Cat. No. GTX64706

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:200
IP	1:20 - 1:50

Not tested in other applications.

Calculated MW 63 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-593 of human ARID3A (NP_005215.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX64706 IHC-P Image

IHC-P analysis of rat lung tissue using GTX64706 DRIL1 antibody.
Dilution : 1:100


GTX64706 IP Image

IP analysis of MCF7 cell lysate using GTX64706 DRIL1 antibody.
Antibody amount : 3 μ g / 150 μ g lysate
Dilution : 1:500


GTX64706 WB Image

WB analysis of various sample lysates using GTX64706 DRIL1 antibody.
Dilution : 1:1000
Loading : 25 μ g per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).