

CD63 antibody [HL5890]

Cat. No. GTX647214

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P, ELISA, Sandwich ELISA
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX646552, Detection: GTX647214

Not tested in other applications.

Product Note This antibody was raised against human CD63 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human CD63. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX647214 ELISA Image

Sandwich ELISA detection of recombinant HEK293 expressed Human CD63 (ECD) protein, human Fc tag (GTX144107-pro) using antibodies as below.

Capture: CD63 antibody [GT368] (GTX646552) (5 µg/mL)

Detection: CD63 antibody [HL5890] (GTX647214) (1 µg/mL)




GTX647214 IHC-P Image

CD63 antibody [HL5890] detects CD63 protein by immunohistochemical analysis.

Sample: Paraffin-embedded human placenta.

CD63 stained by CD63 antibody [HL5890] (GTX647214) diluted at 1:10000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min



GTX647214 IHC-P Image

CD63 antibody [HL5890] detects CD63 protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tonsil.

CD63 stained by CD63 antibody [HL5890] (GTX647214) diluted at 1:10000.

Antigen Retrieval: Tris-EDTA buffer, pH 9.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).