

GPR37L1 antibody [HL6063]

Cat. No. GTX648142

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human GPR37L1 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human GPR37L1. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX648142 IHC-P Image

GPR37L1 antibody [HL6063] detects GPR37L1 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse thalamus.

GPR37L1 stained by GPR37L1 antibody [HL6063] (GTX648142) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX648142 IHC-P Image

GPR37L1 antibody [HL6063] detects GPR37L1 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse cerebellum.

GPR37L1 stained by GPR37L1 antibody [HL6063] (GTX648142) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](#).