

# RBBP9 antibody

## Cat. No. GTX64865

Host	Rabbit	
Clonality	Polyclonal	
Isotype	IgG	
Applications	WB, IHC-P	
Reactivity	Human, Mouse, Rat	

Package 100 μΙ

## Applications

### **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:200

Not tested in other applications.

**Calculated MW** 21 kDa. ( Note )


Properties	
Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-186 of human RBBP9 (NP_006597.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 03 Page 1 of 2


# DATA IMAGES



#### GTX64865 IHC-P Image

IHC-P analysis of human thyroid cancer tissue using GTX64865 RBBP9 antibody.


Dilution: 1:200



### GTX64865 IHC-P Image

IHC-P analysis of rat kidney tissue using GTX64865 RBBP9 antibody.

Dilution: 1:200



### GTX64865 IHC-P Image

IHC-P analysis of human lung cancer tissue using GTX64865 RBBP9 antibody.

Dilution: 1:200



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 03 Page 2 of 2