

Human L-FABP protein

Cat. No. GTX65077-pro

Species Human

Package

25 µg

PRODUCT

Summary Recombinant Human FABP-1 was expressed from E. coli.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Expression System E. coli

Purity ≥98% by SDS-PAGE and HPLC.

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).