

Human Neurotrophin 4 protein (active)

Cat. No. GTX65189-pro

Applications	Functional Assay
Species	Human

Package
10 µg

Applications


Application Note

The ED₅₀ as determined by the dose-dependent induction of choline acetyl transferase activity in rat basal forebrain primary septal cell cultures was found in the range of 20-50 ng/ml.

Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Expression System	E. coli
Purity	??98% by SDS-PAGE and HPLC.
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX65189-pro Image

Purity of recombinant rat Visfatin (2 µg) was analyzed by SDS-PAGE electrophoresis.



For full product information, images and publications, please visit our [website](#).