

DEFB121 antibody

Cat. No. GTX65588

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50 - 1:100

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human DEFB121 (NP_001011878.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX65588 IHC-P Image

IHC-P analysis of human prostate tissue using GTX65588 DEFB121 antibody.
Dilution : 1:100



GTX65588 IHC-P Image

IHC-P analysis of rat testis tissue using GTX65588 DEFB121 antibody.
Blue : DAPI
Dilution : 1:100



GTX65588 IHC-P Image

IHC-P analysis of mouse testis tissue using GTX65588 DEFB121 antibody.
Blue : DAPI
Dilution : 1:100



For full product information, images and publications, please visit our [website](#).