

Human 4-1BBL / CD137L protein, DDDDK tag

Cat. No. GTX65634-pro
Species Human

Package
50 µg

Applications
Product Note Binds to human CD137.

Properties
Form Lyophilized powder

Buffer Reconstitute with 50µl ddH₂O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS (pH7.2).

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence N-terminal DDDDK-Tag; the extracellular domain of human CD137L (a.a. 49-254)


Expression System HEK293 cells

Purity ≥95% by SDS-PAGE.

Endotoxin < 0.1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES

GTX65634-pro Image

"Pull down assay of GTX65634-pro CD137L protein (DDDDK tagged) by hCD137-Fc.

1. hCD137-Fc + hCD137L-DDDDDK-tag
2. hCD137-Fc + hGITRL-DDDDDK-tag
3. hCD137-Fc + hGITRL-His-tag
4. hGITR-Fc + hCD137L-DDDDDK-tag
5. hGITR-Fc + hGITRL-DDDDDK-tag"



For full product information, images and publications, please visit our [website](https://www.genetex.com).