

Human Calreticulin protein, His tag

Cat. No. GTX65637-pro

Species Human

Package

10 µg

Properties

Form Liquid

Buffer 55mM Tris-HCl, 150mM NaCl

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 0.3 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence C-terminal His-Tag; a.a. 18-417

Expression System E. coli

Purity ≥90% by SDS-PAGE.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).