

Human DLK1 protein, DDDDK tag

Cat. No. GTX65645-pro

Species Human

Package
10 µg

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100µl ddH ₂ O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	C-terminal DDDDK-Tag; signal peptide and extracellular domain of human DLK1 (a.a.1-303)
Expression System	HEK293 cells
Purity	≥90% by SDS-PAGE.
Endotoxin	< 0.1 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).