

Human DLK2 protein, human IgG Fc tag

Cat. No. GTX65647-pro

Species Human

Package
10 µg

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence Signal peptide and the extracellular domain of human DLK2 (a.a. 1-306) are fused at the C-terminus to the Fc portion of human IgG.

Expression System HEK293 cells

Purity ≥90% by SDS-PAGE.

Endotoxin < 0.1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).