

# Human Annexin V protein, GST tag

**Cat. No. GTX65648-pro**

**Species** Human

**Package**  
100 µg

## Applications

**Product Note** Binds to phosphatidylserine (PS).

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 100µl ddH<sub>2</sub>O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS (pH7.4).

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** Recombinant human annexin V fused to a GST-tag.

**Expression System** E. coli

**Purity** ≥95% by SDS-PAGE.

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).