

Mouse Periostin protein, DDDDK tag (active)

Cat. No. GTX65650-pro**Applications** Functional Assay**Species** Mouse**Package**

10 µg

Applications

Application Note

Rescues sphere formation from periostin KO primary tumor cells at 50ng/ml.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 100µl ddH₂O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** C-terminal DDDDK-Tag; Mouse periostin (a.a.24-783) (isoform 5)**Expression System** CHO cells**Purity** ??95% by SDS-PAGE.**Endotoxin** < 0.1 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).