

Mouse Periostin protein, DDDDK tag (active)

Cat. No. GTX65650-pro

Applications	Functional Assay
Species	Mouse

Package
10 µg

Applications

Application Note

Rescues sphere formation from periostin KO primary tumor cells at 50ng/ml.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100µl ddH ₂ O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	C-terminal DDDDK-Tag; Mouse periostin (a.a.24-783) (isoform 5)
Expression System	CHO cells
Purity	≥95% by SDS-PAGE.
Endotoxin	< 0.1 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).