

# Mouse IL1 Beta protein

**Cat. No. GTX65651-pro**

**Species** Mouse

**Package**  
10 µg

## Applications

**Product Note** Binds to mouse IL-1R.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 100µl ddH<sub>2</sub>O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** Mouse IL-1 (a.a.118-269)

**Expression System** E. coli

**Purity** ≥95% by SDS-PAGE.

**Endotoxin** < 0.1 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).