

Mouse IL1F8 protein (active)

Cat. No. GTX65654-pro**Applications** Functional Assay**Species** Mouse**Package**

10 µg

Applications

Application Note

Induces secretion of IL-6 in BMDCs at a concentration of 10-100 ng/ml.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 100 µl ddH₂O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS, 0.5% Trehalose.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** a.a. Ser31-Lys183 of Mouse IL-36 beta**Expression System** E. coli**Purity** ≥95% by SDS-PAGE.**Endotoxin** < 0.01 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).