

Mouse IL1F8 protein (active)

Cat. No. GTX65654-pro

Applications	Functional Assay
Species	Mouse

Package
10 µg

Applications

Application Note

Induces secretion of IL-6 in BMDCs at a concentration of 10-100 ng/ml.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100 µl ddH ₂ O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS, 0.5% Trehalose.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	a.a. Ser31-Lys183 of Mouse IL-36 beta
Expression System	E. coli
Purity	??95% by SDS-PAGE.
Endotoxin	< 0.01 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).