

Mouse PEDF protein, DDDDK tag

Cat. No. GTX65656-pro

Species	Mouse
---------	-------

Package	10 µg
---------	-------

Properties

Form	Lyophilized powder
------	--------------------

Buffer	Reconstitute with 100µl ddH ₂ O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
---------	---

Region/Sequence	C-terminal DDDDK-Tag; Mouse PEDF (a.a.20-417)
-----------------	---

Expression System	CHO cells
-------------------	-----------

Purity	≥95% by SDS-PAGE.
--------	-------------------

Endotoxin	< 0.01 EU/µg
-----------	--------------

Conjugation	Unconjugated
-------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---



For full product information, images and publications, please visit our [website](#).