

Toxoplasma gondii Profilin protein

Cat. No. GTX65657-pro

Species Toxoplasma gondii

Package

10 µg

Properties

Form Lyophilized powder

Buffer Reconstitute with 100µl ddH₂O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence Recombinant profilin (Toxoplasma gondii) (a.a. 1-163).

Expression System E. coli

Purity ≥95% by SDS-PAGE.

Endotoxin < 0.01 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).