

Human FNDC5 protein

Cat. No. GTX65658-pro

Species	Human, Mouse, Rat
----------------	-------------------

Package 10 µg

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100µl ddH ₂ O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	FNDC5 (extracellular domain, human (a.a. 32-152) / mouse (a.a. 29-149)
Expression System	E. coli
Purity	≥95% by SDS-PAGE.
Endotoxin	< 0.01 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).