

Mouse CD252 protein, DDDDK tag (active, multimeric)

Cat. No. GTX65660-pro

Applications Functional Assay**Species** Mouse**Package**

10 µg

PRODUCT

Summary

Multimeric Proteins are high activity constructs in which two trimeric TNFSF ligands are linked via the oligomeric collagen domain of ACRP30 [ACRP30headless] and therefore mimic the membrane-bound forms of the proteins.

Applications

Application Note

Activates T cell proliferation.

Product Note

Binds to human and mouse OX40 [CD134].

Properties

Form Lyophilized powder**Buffer** Reconstitute with 100µl ddH₂O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** N-terminal DDDDK-Tag; the extracellular domain of mouse OX40L/CD252 (a.a. 52-199) is fused at the N-terminus to mouse ACRP30headless (a.a. 18-111).**Expression System** HEK293 cells**Purity** ??90% by SDS-PAGE.**Endotoxin** < 0.01 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).