

Human IL1 Receptor 1 protein, human IgG1 Fc tag (active)

Cat. No. GTX65662-pro

Applications Functional Assay

Species Human

Package

50 µg

Applications

Application Note

Blocks IL-1 β activity.

Product Note Binds to human IL-1 β .

Properties

Form Lyophilized powder

Buffer Reconstitute with 50µl ddH₂O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence Human IL-1RI (interleukin-1 receptor type I) (a.a. 21-336) is fused at the C-terminus to the Fc portion of human IgG1.

Expression System HEK293 cells

Purity ??95% by SDS-PAGE.

Endotoxin < 0.1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).