

Human TWEAKR protein, human IgG1 Fc tag

Cat. No. GTX65663-pro

Species Human

Package
50 µg

Applications

Product Note Binds to mouse and human TWEAK.

Properties

Form Lyophilized powder

Buffer Reconstitute with 50 µl ddH₂O to 1.0mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular region of human FN14 (a.a. 28-75) is fused at the C-terminus to the Fc portion of human IgG1.

Expression System HEK293 cells

Purity ≥95% by SDS-PAGE.

Endotoxin < 0.1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).