

Mouse CD24 protein, human IgG1 Fc tag

Cat. No. GTX65664-pro

Species Mouse

Package
50 µg

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 50µl ddH ₂ O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	Mouse CD24 (aa 27-52) is fused at the C-terminus to the Fc portion of human IgG1.
Expression System	HEK293 cells
Purity	≥90% by SDS-PAGE.
Endotoxin	< 0.1 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 11 Page 1 of 1