

Human NAD kinase protein, His tag (active)

Cat. No. GTX65672-pro
Applications Functional Assay

Species Human

Package

50 µg

Applications

Application Note

Biological activity is $\geq 1\text{U/mg}$ protein. One unit is defined as the amount of enzyme that synthesizes $1\mu\text{mol}$ of NADP per min.

Properties

Form Liquid

Buffer 50mM Tris-HCl, 150mM NaCl, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal His-Tag; a.a.1-446

Expression System E. coli

Purity $\geq 95\%$ by SDS-PAGE.

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).