

# Mouse GI24 protein, human IgG1 Fc tag

**Cat. No. GTX65677-pro**

<b>Species</b>	Mouse
----------------	-------

<b>Package</b>	50 µg
----------------	-------

**Properties**

<b>Form</b>	Lyophilized powder
-------------	--------------------

<b>Buffer</b>	Reconstitute with 50µl ddH <sub>2</sub> O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.
---------------	---

<b>Preservative</b>	No preservatives
---------------------	------------------

<b>Storage</b>	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

<b>Region/Sequence</b>	The extracellular domain of mouse VISTA (a.a. 33-191) is fused at the C-terminus to the Fc portion of human IgG1.
------------------------	---

<b>Expression System</b>	HEK293 cells
--------------------------	--------------

<b>Purity</b>	≥95% by SDS-PAGE.
---------------	-------------------

<b>Endotoxin</b>	< 0.01 EU/µg
------------------	--------------

<b>Conjugation</b>	Unconjugated
--------------------	--------------

<b>Note</b>	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---



For full product information, images and publications, please visit our [website](#).