

# Human Jagged 2 protein, human IgG Fc tag

**Cat. No. GTX65678-pro**

**Species** Human

**Package**

10 µg

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 100 µl ddH<sub>2</sub>O to 0.1mg/ml. Lyophilized from 0.2µm-filtered PBS, 0.02 % CHAPS.

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** Human Jagged-2 (a.a. 1-1054) is fused at the C-terminus to the Fc portion of human IgG.

**Expression System** HEK293 cells

**Purity** ≥90% by SDS-PAGE.

**Endotoxin** < 0.1 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).