

Human DR6 protein, human IgG1 Fc tag

Cat. No. GTX65679-pro

Species	Human
----------------	-------

Package

50 µg

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Reconstitute with 50µl ddH ₂ O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.
---------------	---

Preservative	No preservatives
---------------------	------------------

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

Region/Sequence	The Fc portion of human IgG1 is fused to the C-terminus of human DR6 (a.a. 42-351).
------------------------	---

Expression System	HEK293 cells
--------------------------	--------------

Purity	≥95% by SDS-PAGE.
---------------	-------------------

Endotoxin	< 0.1 EU/µg
------------------	-------------

Conjugation	Unconjugated
--------------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---



For full product information, images and publications, please visit our [website](#).