

Human CD276 protein (active)

Cat. No. GTX65682-pro

Applications	Functional Assay
Species	Human

Package
50 µg

Applications

Application Note

Measured by its ability to inhibit anti-CD3-induced proliferation of stimulated human T cells.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100µl ddH ₂ O to 0.5mg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	The extracellular domain of human CD276 [B7-H3] (a.a. 29-245)
Expression System	HEK293 cells
Purity	??95% by SDS-PAGE.
Endotoxin	< 0.01 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).