

Human PD-L1 protein (active)

Cat. No. GTX65685-pro

Applications Functional Assay

Species Human

Package

50 µg

Applications

Application Note

Measured by its ability to inhibit anti-CD3 antibody induced IL-2 secretion in human T lymphocytes.

Properties

Form Lyophilized powder

Buffer Reconstitute with 100µl ddH₂O to 0.5mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular domain of human CD274 [B7-H1] (a.a.19-238)

Expression System HEK293 cells

Purity ≥95% by SDS-PAGE.

Endotoxin < 0.01 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).