

# Human ICOS Ligand protein, His tag

**Cat. No. GTX65686-pro**

**Applications** Functional Assay

**Species** Human

**Package**

50 µg

## Applications

### Application Note

Measured by its binding ability in a functional ELISA.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 100µl ddH<sub>2</sub>O to 0.5mg/ml. Lyophilized from 0.2µm-filtered PBS.

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** C-terminal His-Tag; a.a. 19-256

**Expression System** HEK293 cells

**Purity** ??95% by SDS-PAGE.

**Endotoxin** < 0.01 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).