

Human Gl24 protein, His tag

Cat. No. GTX65689-pro

Species Human

Package

50 µg

Properties

Form Lyophilized powder

Buffer Reconstitute with 100µl ddH₂O to 0.5mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence C-terminal His-Tag; a.a. 33-194

Expression System HEK293 cells

Purity ≥95% by SDS-PAGE.

Endotoxin < 0.01 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).