

Human CD40L / CD154 protein, human IgG1 Fc tag

Cat. No. GTX65706-pro**Species** Human**Package**

50 µg

Properties

Form Lyophilized powder**Buffer** Reconstitute with 50 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of human CD40L [CD154] (a.a. 113-261) is fused to the N-terminus of the Fc region of human IgG1.**Expression System** CHO cells**Purity** ≥98% by SDS-PAGE.**Endotoxin** < 0.06 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).