

# Human Galectin 9 protein, human IgG1 Fc tag

**Cat. No. GTX65711-pro**

**Applications** Functional Assay

**Species** Human

**Package**

50 µg

## Applications

### Application Note

Measured by its binding ability in a functional ELISA.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 50 µl ddH<sub>2</sub>O. Lyophilized from 0.2µm-filtered PBS.

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** The extracellular domain of human Galectin-9 (a.a. 2-323) is fused to the N-terminus of the Fc region of human IgG1.

**Expression System** HEK293 cells

**Purity** ??95% by SDS-PAGE.

**Endotoxin** < 5 EU/mg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).